Intelligent
evaporation

Betriebsanleitung
Operating instructions
Notice d'instructions

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Translation of the original instructions
Page Page 80 – 148

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Contents

Introduction

About this document	80
Typographic conventions	80
Copyright protection	80

General information

Notes on the product	81
Directives applied, product certification	81
California Residents	81
Customs declaration	81
Residual risk	81
Intended use	81
Reasonably foreseeable misuse	82
Transportation	82
Storage	82
Acclimatization	82
Permissible ambient conditions	82

Safety

General safety instructions	83
Electrical safety	83
Data security	83
Operational safety	83
Work safety	84
Personal protective equipment (PPE)	84
Environmental protection	84
Biohazard	84
Special hygiene measures for the use of laboratory equipment in food, cosmetics and pharmaceutical production	85
General measures	85
Device-specific measures	85
Other regulations	86

Mounting the system

Mounting	87
----------------	----

Device description

Mechanical design	88
Device connections	89
Control-Box	90
Interfaces.....	90
Control panel	91
Structure of the user interface.....	92
Control and display elements.....	94
Buttons and symbols	95
On-screen keyboard	96

Commissioning

Commissioning.....	97
Connecting the control panel	97
Power supply.....	97
Switching the device on/off	98
Transportation lock.....	98
Evaporation flask	99
Mount/remove evaporation flasks.....	99
Evaporation flask inclination	100
Evaporation flask immersion depth.....	101
Flask lift	102
Receiving flask	103
Heating bath.....	104
Heating bath	104
Filling the heating bath	105
Emptying the heating bath	105
Heating bath overheat protection	106
Residual heat indicator.....	106
Device configuration	107
System time, timer, stopwatch.....	107
Main menu.....	109
Applications	111
Settings.....	112
Selecting the user language	118
Selecting formats.....	119
Data Recording	120
Activating/deactivating the data logging function.....	120

Connectivity123
System check124
Cleaning the PTFE seal125

Operation

Operating modes126
Base mode127
Advanced mode.....128
 Dynamic AUTOaccurate (DAA) mode..... 128
 Ramp profiles mode 131
 Favorites mode 137

Troubleshooting

Troubleshooting.....141

Appendix

Technical specifications.....143
Scope of delivery144
Accessories.....144
Device service.....145
 General cleaning instructions 145
 Repairs 145
 Maintenance..... 146
 Dismantling..... 146
 Disposal..... 146
Contact information Heidolph international.....147
Warranty Statement147
Certificate of decontamination148

About this document

These operating instructions describe all the functions and operation of the Hei-VAP Expert, Hei-VAP Expert Control, Hei-VAP Ultimate and Hei-VAP Ultimate Control type rotary evaporators. The operating instructions are an integral part of the described device!

Typographic conventions

Standardized symbols, signal words and highlighting are used in this document to warn of hazards, and to identify important information and special text contents.

Symbol	Signal word / explanation
	<p>Warning symbols in combination with a signal word indicate dangers:</p> <p>DANGER Indicates an immediate dangerous situation which, if not avoided, will result in death or serious injury.</p> <p>WARNING Indicates a potential danger which, if not avoided, may result in serious injury.</p> <p>CAUTION Indicates a potential hazard which, if not avoided, can result in damage to property and minor to moderate injuries.</p>
	<p>Mandatory signs indicate important and useful information on handling a product.</p> <p>This information is used to ensure operational safety and to maintain the value of the product.</p>
[GUI]	<p>Parameter designations, display texts, and device labels are typographically highlighted in the text and in tables to facilitate their assignment on the device.</p>
→	<p>The arrow indicates specific instructions to be followed to ensure operational safety when handling the product.</p>

Copyright protection

This document is protected by copyright and is intended for use by the purchaser of the product only.

No transfer to third parties, reproduction in any form, including excerpts, and by any means, and utilization and/or disclosure of the contents is permitted without the prior written consent of Heidolph Instruments GmbH & Co. KG. Any violation is subject to compensation for damage.

Notes on the product

Directives applied, product certification

	<p>CE marking</p> <p>The device meets all requirements of the following directives:</p> <ul style="list-style-type: none"> • European Machinery Directive, 2006/42/EC • EMC Directive, 2004/108/EC
	<p>NRTL Certification</p> <p>The device has been tested in accordance with the following standards:</p> <ul style="list-style-type: none"> • UL 61010-1 :2012/R:2016-04 CAN/CSA-C22.2 No. 61010-1:2012/U2:2016-04 • UL 61010-2-010:2015 CAN/CSA-C22.2 NO. 61010-2-010:2015

California Residents

Important information for California residents regarding Prop 65. Please visit www.P65Warnings.ca.gov for more information.

Customs declaration

Rotary evaporators, as combustion and cleaning equipment, may be subject to notification to a competent customs authority in the country of destination. The evaluation of a customer obligation to notify and, if applicable, the notification to a competent customs authority in the country of destination is solely the responsibility of the user!

Residual risk

The device has been designed and manufactured in accordance with the state-of-the-art standards at the time of development and the recognized safety regulations. During mounting and use, as well as during maintenance, repair and cleaning work, there are nevertheless certain residual risks associated with the described device.

These are identified and described at the appropriate points in this document.

Intended use

The described device was designed by the manufacturer for the separation of substances, powder drying, concentration, crystallization of substances, and recycling of solvents under vacuum.

Due to its design, the device in its delivery condition may only be used in analytical processes or in laboratory-like conditions in the food, cosmetics, and pharmaceutical industries as well as other comparable industries that manufacture products intended for consumption by humans or animals, or for use on humans or animals.

Any other use of this device is not considered as intended!

Reasonably foreseeable misuse

Additional measures may be necessary, and/or specific directives and safety regulations may have to be observed for use under conditions or for purposes deviating from the intended use (see, for example, section "Special hygiene measures for the use of laboratory equipment in food, cosmetics and pharmaceutical production" on page 85). Corresponding requirements must be evaluated and implemented by the operator in each individual case.

Compliance with and implementation of all relevant directives and safety measures for the respective field of application is the sole responsibility of the operator.

All risks resulting from improper use are borne solely by the operator.

The device may only be operated by authorized and instructed personnel. Training and qualification of the operating personnel as well as ensuring that the device is handled responsibly are the sole responsibility of the operator!

Transportation

During transportation, avoid severe shocks and mechanical stresses that can cause damage to the device.

Keep the original packaging in a dry and protected place for later use.

Storage

Always store the device in its original packaging. To protect against damage and excessive material aging, store the device in an environment that is as dry, temperature-stable and dust-free as possible.

Recommended storage conditions:

- 5 °C – 31 °C up to 80 % rel. humidity
- 32 °C – 40 °C up to 50 % rel. humidity (decreasing linearly)

Acclimatization

After each transportation and after storage under critical climatic conditions (e.g. high temperature difference between inside and outside), allow the device to acclimatize at room temperature for at least two hours to prevent possible damage due to condensation before commissioning it in the place of use. If necessary, extend the acclimatization phase if the temperature differences are very high.

Make all supply connections (power supply, tubing) only after the device has been acclimatized!

Permissible ambient conditions

The device is designed for indoor use only. Permissible ambient conditions for operation:

- 5 °C – 31 °C at up to 80 % rel. humidity
- 32 °C – 40 °C at up to 50 % rel. humidity (decreasing linearly)
- Installation altitude: 2,000 m asl

When used in corrosive atmospheres, the service life of the device may be reduced depending on the concentration, duration and frequency of exposure.

The device **IS NOT** suitable for outdoor use!

The device **IS NOT** suitable for use in potentially explosive areas!

General safety instructions

- Before commissioning and using the device, familiarize yourself with all the safety regulations and occupational safety guidelines applicable at the place of use and observe them at all times.
- Only operate the device if it is in faultless technical condition. In particular, ensure that there is no visible damage on the device itself, the power supply cord and, where applicable, on connected devices or the supply connections.
- If there is missing or misleading information on the device or regarding occupational safety, contact the responsible safety specialist or our technical service.
- Only use the device in accordance with the regulations on intended use ("Intended use" on page 81).

Electrical safety

- Before connecting the device to the power supply, ensure that the voltage indicated on the rating plate matches the specifications of the local power utility company.
- Ensure that the power supply circuit provided is protected by means of a residual-current device (RCD).
- Only use the three-pole power supply cord supplied with the device.
- The device must only be supplied with power from a properly grounded mains socket-outlet.
- Have repairs and/or maintenance work on the device carried out exclusively by an authorized electrician or by the technical service department of Heidolph Instruments.
- Always disconnect the device from the power supply system before carrying out any maintenance, cleaning or repair work.

Data security

- The user is responsible for ensuring data security when transferring data between the described device and other devices.
- Only use secure networks for the data transfer and avoid use of critical infrastructure.
- Only use high-quality shielded data cables for the data transfer.
- For data transfer via a USB B connection, an industrial standard USB hub should be preferably used to ensure the most stable connection possible.

Operational safety

- Operate the device under a closed ventilated fume hood when working with potentially hazardous substances (see EN 14175 and DIN 12924).
- Do not make any unauthorized changes or modifications to the device!
- Only use genuine spare parts and accessories, or those expressly approved by the manufacturer!
- Rectify malfunctions or faults on the device immediately.
- Switch off and disconnect the device from the power supply, preventing reconnection, if it is not possible to eliminate the malfunction or rectify the fault immediately.
- Observe all other applicable regulations such as laboratory and workplace guidelines, recognized safety technology rules and special local regulations.

Work safety

- Always use the prescribed personal protective equipment (PPE) such as protective clothing, safety goggles, protective gloves, safety shoes, etc.
- Persons with pacemakers or other susceptible implants must maintain a minimum distance of 10 cm from the device to prevent possible interference.
- Do not operate any other devices in the immediate vicinity of the device,
 - which generate emission or radiation sources in the frequency range 3×10^{11} Hz to 3×10^{15} Hz (in the optical spectral range wavelengths from 1,000 μm to 0,1 μm),
 - which generate ultrasonic or ionizing waves.
- Do not process any substances that could release energy in an uncontrolled manner (e.g. self-ignition).
- Do not process substances in which the energy input through mixing poses a danger.
- Do not work with an open flame in the immediate vicinity of the device (risk of explosion).
- Wipe off any fluid that may have spilled on the device immediately.
- Always switch off the device when it is not in use.

Personal protective equipment (PPE)

The operator must determine and provide the necessary PPE, depending on the respective application and the media and chemicals used.

The corresponding instruction of the personnel is solely within the operator's responsibility.

Environmental protection

When processing environmentally hazardous substances, take appropriate measures to avoid risks to the environment.

The evaluation of corresponding measures such as the marking of a hazardous area, their implementation, and the training of the relevant personnel is the sole responsibility of the operator!

Biohazard

When processing biohazardous substances, take appropriate measures to prevent hazards to persons and the environment, including:

- Instruction of the personnel regarding the necessary safety measures.
- Provision of personal protective equipment (PPE) and instruction of the personnel in its use.
- Marking the device with a biohazard warning symbol.

The evaluation of corresponding measures such as the marking of a hazardous area, their implementation, and the training of the relevant personnel is the sole responsibility of the operator!

Special hygiene measures for the use of laboratory equipment in food, cosmetics and pharmaceutical production

When laboratory equipment is used in the production processes of the food, cosmetics or pharmaceutical industry, special hygiene measures must be taken by the user to avoid sample contamination and to minimize any risk to humans and the environment as far as possible.

Please observe the following manufacturer's recommendations:

General measures

- Ensure a clean working and storage environment when handling substances and materials.
- Train all employees in the field of occupational hygiene, document all training measures and check the implementation of all required hygiene measures regularly during operation.
- Use a hygiene control concept such as HACCP (Hazard Analysis and Critical Control Points). The HACCP comprises the following criteria:
 - Hazard analysis
 - Identification of critical control points
 - Definition of critical limit values
 - Establishment of a system for monitoring and controlling critical control points (CCP)
 - Corrective actions for uncontrollable CCP
 - Establishment of a system to verify the implementation of all HACCP measures
 - Establishment of a system for documenting all associated procedures and records

The evaluation of the applicability of the mentioned rules and regulations is within the sole responsibility of the operator!

Device-specific measures

- Regularly clean components that come into contact with the product, such as flasks, seals, tubes, etc. in the autoclave (if available or possible) or chemically (e.g. with ethanol) to sterilize all surfaces.
- Make sure that even products that are intended for single use only are of sufficient purity.
- Do not use open containers.
- Avoid contamination caused by careless handling of contaminated vessels, apparatus or aids.

Contact information

For further information, please contact our after sales service at any time.

Phone: +49-9122-9920-0

Email: sales@heidolph.de

Other regulations

In addition to the notes and instructions in this document, observe all other applicable regulations such as laboratory and workplace guidelines, hazardous substances ordinances, recognized rules of safety engineering and occupational medicine as well as particular local regulations!

Noncompliance will invalidate any warranty claims against Heidolph Instruments.

The operator is solely liable for all damage resulting from unauthorized changes or modifications to the device, from the use of unauthorized or non-genuine spare parts and accessories, or from disregarding the safety instructions and hazard warnings or the manufacturer's instructions!

Mounting

If required, Heidolph offers a professional system and/or device mounting service, including all supplied components and including commissioning.

The associated mounting instructions describe all steps for proper system and/or device mounting in detail!

CAUTION: Damage to the device

Due to improper mounting and/or installation of the system/device, as well as unauthorized changes to the mounting, there is a risk of direct and indirect damage to property!

- Do not make any unauthorized changes to the system/device mounting.
 - Do not make any unauthorized changes to the tubing and power supply connections.
 - If necessary (system mounting and installation carried out by the operator or third parties, necessary changes to the existing mounting), contact the manufacturer's technical service (see section "Contact information Heidolph international" on page 147).
-

Ensuring proper system and/or device mounting, including all accessory parts, is the sole responsibility of the operator!

Noncompliance will invalidate any warranty against Heidolph Instruments.

The operator is solely liable for all damage resulting from unauthorized changes or modifications to the device, from the use of unapproved or non-genuine spare parts and accessories, or from disregarding the safety instructions and hazard warnings or the manufacturer's instructions!

Mechanical design

The Hei-VAP Expert and Hei-VAP Ultimate variants have an identical design, but are equipped with different control panels. The following illustration shows an example of the Hei-VAP Ultimate model with glassware G3:


Control panel
Hei-VAP Expert


Device connections

Hei-VAP Expert Control and Hei-VAP Ultimate Control with vacuum valve for use with diaphragm vacuum pump (Hei-VAP Expert and Hei-VAP Ultimate, each without Control Box!):

Hei-VAP Expert Control and Hei-VAP Ultimate Control for use with speed-controlled vacuum pump (Hei-VAP Expert and Hei-VAP Ultimate, each without Control Box!):

Control-Box

The Control Box is an integral part of the delivery of the Hei-VAP Expert Control and Hei-VAP Ultimate Control variants.

The Control Box is used to control the peripheral devices.

For detailed information on connecting the associated system components, see the installation instructions delivered with each unit!

Interfaces

The device interfaces are located on the side of the device display:

- 1 × USB 2.0, max. 500 mA (USB type interface: software and firmware updates, data backup)
- 1 × Micro SD (interface for service purposes)

CAUTION: Damage to the device

Make sure that the device interfaces on the bottom of the display are always properly covered!

Penetrating fluid can cause interface malfunctions and short-circuit.

Control panel

The control panels of the Hei-VAP Ultimate and Hei-VAP Expert variants differ in terms of the dimensions and arrangement of the controls. The functionality of the controls described is identical for all variants.

The control elements of the central display are touch-sensitive (touch panel) and can be operated by finger or with a suitable stylus.

CAUTION: Damage to the control panel

Tip or sharp objects can damage the surface of the display on the control panel!


→ To operate the device, touch the display only with your fingertips or a special stylus with a soft-touch tip.

Structure of the user interface

All device functions are controlled via the integrated control panel and the graphical user interface.


The displays of the variants Hei-VAP Ultimate and Hei-VAP Expert differ in terms of the dimensions and arrangement of the controls, the functionality of the described elements is identical:

For clarity, all controls in the following sections are described using the example of the Hei-VAP Ultimate variant.

Where applicable, deviations of the Hei-VAP Ultimate variant are clearly identified and described.

After switching on the device and initializing the device control, the start window of the last active operating mode appears on the display (in the example, DAA automatic mode):

From the start window, you can access the main menu, select an operating mode directly, and select the time display, see also "Control and display elements" on page 94.

Control and display elements

Control element	Function
Main menu function button	<p>Touch this function button to open the main menu with the following options:</p> <ul style="list-style-type: none"> ▪ Applications → operating mode ▪ Settings → system parameters ▪ Languages → menu language ▪ Formats → system configuration and display formats ▪ Data logging → data recording ▪ System check → complete system check ▪ Error list → open error list ▪ Info → open system info window
Function button Operating mode	<p>Touch this function button to activate the required operating mode:</p> <ul style="list-style-type: none"> ▪ Automatic mode ▪ Ramp profile mode ▪ Favorites mode
Function button System time/Timer	<p>Touch this function button to toggle between the following functions:</p> <ul style="list-style-type: none"> ▪ System time display ▪ Timer ▪ Stopwatch <p>See the following section "System time, timer, stopwatch" on page 107 for detailed information.</p>
Start process button	Touch this button to start the selected process.
Warm-up button	Touch this button to warm up the heating bath.
Stop process button	Touch this button to stop the running process.
Set value/actual value indicator process parameters	<p>The following process parameters are displayed:</p> <ul style="list-style-type: none"> ▪ Boiling temperature/AUTOaccurate temperature ▪ Coolant temperature ▪ Vacuum pressure ▪ Rotation speed ▪ Heating bath temperature
Process visualization main display	The current process step is visualized in the main display.

Buttons and symbols


The following table provides an overview of the most important HMI buttons and icons on the control panel and in the menu navigation:

Symbol	Function
	Menu button
	Home button
	Automatic mode DAA button
	Ramp profile mode button
	Favorites mode button
	Start / start all button
	Stop / stop all button
	Confirm/OK button
	Cancel/Discard/NOK button
	Delete button
	Screen lock icon

On-screen keyboard

The numeric on-screen keyboard appears automatically when you touch a numeric entry field within a mask.

The alphanumeric on-screen keyboard appears automatically when you touch an input field for alphanumeric values within a mask.

- Use the keyboard to enter the desired designation or value within the permitted range of values.
- Touch the **Cancel** button to return to the parent dialog without making any changes.
- Touch the Delete button to cancel the last entry.
- Confirm your entry with **ok**. The on-screen keyboard closes.

Commissioning

Follow the steps below to start up the Hei-VAP Expert or Hei-VAP Ultimate rotary evaporator.

Connecting the control panel

Before commissioning the device, connect the cable of the removable control panel to the cable connector on the base unit:

Power supply

DANGER: Electric Shock

- Only use the power supply cord provided with the device!
- Before connecting the device to the power supply, ensure that the voltage indicated on the rating plate matches the specifications of the local power utility company.
- Ensure that the power supply circuit provided is protected by means of a residual-current device (RCD).
- Switch off the device immediately and secure it against being switched on again or do not put the device into operation if the power supply cord or the associated connection points show visible damage or corrosion. In such a case, contact a qualified and authorized electrician!

- Connect the power supply cord to the appliance inlet on the back of the base unit.
- Make sure that the device is switched off (main switch in position off, see section "Mechanical design" on page 88).
- Connect the power supply cord to a properly secured mains socket-outlet.

Switching the device on/off

CAUTION: Risk of property damage, loss of production

Before each operation of the device, check the sealing caps and connecting elements of all glass components for correct and secure positioning!

At loose/leaky connections, ambient air is sucked in during the vacuum build-up and the required vacuum pressure cannot be reached in a stable manner!

In case of excess pressure, there is a risk that fluid will escape from the system! The device and surrounding areas can be severely contaminated by leaking media!


→ Tighten loose caps/fasteners by hand or correct the seating of the retaining clips.

To switch on and off, use the main switch located at the front side of the device, see section "Mechanical design" on page 88.

Transportation lock

All variants with motor driven lift are provided with a mechanical transportation lock to prevent damage during the transport of the device. This transportation lock must be removed before commissioning the device! Keep the transportation lock with the original packaging for future use.

The transportation lock consists of a profiled plate, which is fixed to the housing with three screws M5×8. The profiled plate is marked with a red flag:

To remove the transportation lock, follow these steps:

- Make sure that the heating bath and the control panel are connected correctly (see sections "Connecting the control panel" on page 97 and "Heating bath" on page 104).
- Switch on the device (see section "Follow the steps below to start up the Hei-VAP Expert or Hei-VAP Ultimate rotary evaporator." on page 97).
- Using the supplied hex wrench, loosen the three fixing screws and remove the profiled plate.
- Move the flask lift to the upper end position (see section "Flask lift" on page 102).

- Turn off the device.
 - Switching off the device with the flask lift in the upper end position, the additional electronic safety function will be disabled.
- Switch on the device again.
- The device is then ready for use.

Evaporation flask

CAUTION: Risk of scalding, burns and/or damage to property

In case of improper handling, the evaporation flask may fall to the ground.

There is a risk of glass breakage.

Residues remaining in the flask can run out and cause burns, scalding and/or damage to property, as well as production downtime.

- Refer to the instructions in this section when mounting/removing an evaporation flask.
- Use only the evaporating flasks approved by the manufacturer.

WARNING: Risk of implosion, risk of poisoning

The smallest damage to the evaporation flask can lead to glass breakage when pressurized. As a result, there is a risk of injury due to glass splinters and the risk of uncontrolled release of distillation material.

- Before each operation of the device, check the evaporation flask for visible damage.

Mount/remove evaporation flasks

The evaporation flask is fixed to the vapor tube of the rotary evaporator using a special retort clamp (EASY CLIP):

EASY CLIP open, flask unclamped

EASY CLIP closed, flask clamped

To mount and remove the evaporation flask, follow these steps:

- Stop all running processes and the rotation of the flask.
- Move the flask lift to the upper end position (see section "Flask lift" on page 102).
- Loosen the screw connection of the EASY CLIP at the vapor tube and push the clamping bracket outward (see figure above).
- When mounting the flask, make sure that the glass cones and sockets of the flask and the vapor tube do have full contact and form a straight line (clean the glass cones/sockets before mounting the flask!).
- Clamp the flask by closing the bracket of the EASY CLIP until it arrests with an audible click.
- Screw the EASY CLIP back onto the vapor tube.
- To immerse the flask in the heating bath, move the flask lift to the lower end position (see section "Flask lift" on page 102). Take into consideration the flask size and the filling quantity of the heating bath (see section)! "Filling the heating bath" on page 105

CAUTION: Risk of scalding

The heating bath is designed for a maximum operating temperature of 210°C (using silicone oil). From a liquid temperature of 50 °C, there is a risk of injury in the event of contact.

- Allow heated liquid to cool below 50 °C before removing the flask.

or

- Use appropriate personal protective equipment (heat-resistant gloves, eye protection, safety clothing) to remove the flask.

Evaporation flask inclination

To adjust the inclination angle of the evaporation flask, do the following:

- Stop all running processes and the rotation of the flask.
- Move the flask lift to the upper end position (see section "Flask lift" on page 102).
- Loosen the locking on the lift base and hold it against the spring force.
- Swing the condenser until the evaporation flask is at the desired angle. Make sure that the flask can turn freely!

- Release the locking and if necessary, slightly move the unit condenser-evaporator flask until the drive unit clicks into place.

Evaporation flask immersion depth

To adjust the immersion depth of the evaporation flask, do the following:

- Stop all running processes and the rotation of the flask.
- Move the flask lift with the evaporation flask mounted to the desired height (see section "Flask lift" on page 102).
- Loosen the fixing screw and push the stop adjustment bar down to the lift stop (see the following figure).
- Tighten the fixing screw again.
- If necessary, position the base plate with the heating bath according to the flask size.

The lift can only be moved downwards up to the set immersion depth. Before the evaporation flask is immersed, check that the fixing screw is firmly seated!

When adjusting the immersion depth, also observe the instructions in section "Filling the heating bath" on page 105 to avoid the heating bath fluid overflowing!

Flask lift

Devices with manual lift

Devices with a manual lift are equipped with a lever for moving the lift up and down:

- Push and move the lever to the left to move the lift up.
- Push and move the lever to the right to move the lift down.

Devices with motor driven lift

For devices with motor driven lift, the lift is moved using the two lift buttons on the device display:

- Touch the UP button to move the lift upwards.
- Touch the DOWN button to move the lift downwards.
- The movement stops as soon as you release the button.

Receiving flask

The receiving flask is fixed to the condenser of the rotary evaporator with a retort clamp:

To mount and remove the receiving flask, follow these steps:

- Stop all running processes and the rotation of the flask.
- Move the flask lift to the upper end position (see section "Flask lift" on page 102).
- While holding the receiving flask, loosen (turn the knurled screw counter-clockwise) and pull off the ground clamp (if applicable).
- When mounting the receiving flask, make sure that the glass cones and sockets of the flask and the condenser do have full contact and form a straight line (clean the glass cones/sockets before mounting the flask!).
- Mount the ground clamp over the joint and tighten the ground clamp by hand (turn the thumbscrew clockwise).

Heating bath

Observe the instructions in this section when handling the heating bath and the bath fluid.

Heating bath

Always grip and hold the heating bath by the ergonomic safety handles on the side!

The heating bath is placed on the flexible base plate of the base unit and connected to the unit by means of a 7-pin connecting cable. The recess provided on the bottom side ensures that the heating bath can only be placed in the right position on the base plate.

7 pole connection cable
heating bath - base unit

To place the heating bath on the base unit, follow these steps:

- Stop all running processes and the rotation of the flask (if applicable).
- Move the flask lift to the upper end position (see section "Flask lift" on page 102).
- Place the heating bath on the base plate of the basic unit: The base plate is equipped with a guide suitable for the bottom-side recess of the heating bath.
- Make sure that the heating bath is securely fixed to the base plate and in position.
- Connect the supplied 7-pin connection cable between the socket of the heating bath and the base unit.
- Before removing the heating bath, disconnect the 7-pin connecting cable between the heating bath socket and the base unit!

Filling the heating bath

CAUTION: Risk of property damage, loss of production

If the filling level is too high, there is a risk that the heating bath fluid will overflow when the flask is immersed.

The device and surrounding areas can be severely contaminated by escaping heating bath fluid!

- Take the imprinted MIN/MAX markings on the inside of the heating bath as a reference.

- Ensure that the evaporation flask is correctly mounted ("Mount/remove evaporation flasks" on page 99).
- Move the flask lift to the lower end position (see section "Flask lift" on page 102).
- With the flask in the lower end position, fill the heating bath until the fluid reaches a level between the MIN/MAX marks on the inside.

Emptying the heating bath

CAUTION: Risk of scalding

The heating bath is designed for a maximum operating temperature of 210°C (using silicone oil). From a liquid temperature of 50 °C, there is a risk of injury in the event of contact.

- Always grip the heating bath by the ergonomic safety handles.
- Before emptying the heating bath, allow heated liquid to cool down to below 50 °C.

or

- Use suitable personal protective equipment (heat-resistant gloves, eye protection, safety clothing) to empty the heating bath.

- Move the flask lift to the upper end position (see section "Flask lift" on page 102).
- Disconnect the 7-pin connecting cable between the heating bath and the base unit.
- Remove the heating bath from the base plate of the base unit.
- Empty and clean the heating bath.

When changing the bath fluid, the heating bath must be thoroughly cleaned and dried.

In particular, when changing from water to oil, make sure that the heating bath temperature is only slowly and gradually augmented during restart (recommendation: 10 °C), so that remaining water in the system can evaporate until the process temperature is reached.

Heating bath overheat protection

The device is equipped with an autonomous overheat protection function that prevents the heating bath from overheating. When reaching the maximum heating bath temperature, the heating function is switched off by a mechanical switch

To reset the safety function, follow these steps:

- Remove the heating bath from the base unit.
- Empty and let cool down the heating bath.
- Press the reset button on the bottom of the heating bath with a suitable tool (blunt tip, insulated or non-conductive material).

- Put the heating bath back on the base unit.


Residual heat indicator

As long as the heating bath temperature is above 50° C, the LED ring of the temperature/vacuum control knob flashes after the end of the process. If the temperature/vacuum control knob is assigned the vacuum control function, a temperature warning message will appear on the display instead.

The residual heat indicator goes out when the device is switched off!

Device configuration


The basic settings for operation are made in the main menu. After switching on the device and initializing the software, the start page of the last active operating mode appears (in the example: Advanced mode, DAA mode). To open the main menu, touch the **Main menu** function button:

System time, timer, stopwatch

Touch the **System time/Timer** function button in the home screen header to open the **time/timer/stopwatch** window. Depending on the mode selected, the following information is displayed:

- Current time (12- or 24-hour display, depending on the selection in the **Settings** menu).
- Stopwatch: Display of the actual elapsed time from the start of the process.
- Timer: Display of the time remaining until the scheduled end of the process.

This setting is applied in all operating modes and can be changed at any time!

Setting the date and time

- Switch to **Time** mode.
- Touch and hold the **System time/Timer** function button for about two seconds to open the system time setting dialog box:


- Touch the **Date** function button and use the keyboard to enter the current date.
- Touch the **Time Zone** function button and select your time zone.
- Touch the **Time** function button and use the keyboard to enter the current time.

Programming the timer

- Switch to **Timer** mode.
- Touch the **System time/Timer** function button for approx. two seconds to open the on-screen keyboard (see section "On-screen keyboard" on page 96) and to program the timer.

Main menu


The following options are available from the main menu:

Button	Function/Meaning
	<p>Home button. Touch this function button to return to the previously opened view (start window).</p>
Operating modes	<p>Touch this function button to open the Applications menu and select an operating mode, see section "Applications" on page 111.</p>
Settings	<p>Touch this function button to open the Settings menu. In this menu, various basic settings can be adjusted and the device can be reset to the factory settings, see section "Settings" on page 112.</p>
Languages	<p>Touch this function button to open the Languages menu and load a required user language, see section "Selecting the user language" on page 118.</p>
Formats	<p>Touch this function button to open the Formats menu. Here you can select the desired format for the display of numerical values, time, date and temperature. At the same time, the connected rotary evaporator / glassware is assigned and the operation with or without residue drainage is selected, see section "Selecting formats" on page 119.</p>
Data logging	<p>Touch this function button to start the Data logging function, see section "Data Recording" on page 120.</p>
Hei-CONTROL Pro	<p>Touch this function button to open the screen form for definition of the network parameters, see section "Connectivity" on page 123.</p>
System check	<p>Touch this function button to perform a full system check, see section "System check" on page 124.</p>
Error list	<p>Touch this function button to open the list of all error messages.</p>
Info	<p>Touch this function button to open the information window of the system, including firmware and hardware data .</p>

Applications

In the **Applications** menu (**Main menu** → **Applications**), you can switch between manual mode (**Basic**) and automatic operation with ramp profiles or favorites (**Advanced**) (see section "Operating modes" on page 126 for a detailed description of the operating modes):


To activate the desired operating mode, proceed as follows:

- Complete all running processes before changing the operating mode.
- Touch the entry of the desired operating mode. The background color of the selected entry changes to orange, as shown.
- Press **ok** (green checkmark) to confirm the selection and activate the selected mode.
 - The start window of the activated mode appears on the display.
- Press **discard** (red X) to discard the selection and to return to the previous view.

Settings

In the **Settings** menu (**Main menu** → **Settings**), various device settings can be adjusted:

The menu includes several screens! Touch the display and swipe up or down to scroll through the menu. The operating functions are identical in all screens:

- Touch the **Cancel** button to return to the previous view or the start window without making any changes.
- Confirm your entries with **ok** to accept changes and to return to the previous view or the start window of the operating mode.

The various options are described in detail in the following sections:

Menu item **Rotary knob right**:

Use this menu item to define the function of the temperature/vacuum control knob on the control panel.

- Touch the function button to switch between the two options:
 - Selection **Vacuum**: Knob is used to regulate the vacuum pressure.
 - Selection **Heating bath**: Knob is used to regulate the heating bath temperature.

Menu item **Hysteresis** Δ

In this menu, define the hysteresis for the opening pressure of the vacuum valve.

- Touch the function button and choose between automatic and manual hysteresis mode:
 - Activating the **Auto** mode, predefined hysteresis values are applied for different pressure ranges.
 - Select **Manual hysteresis** to define a uniform hysteresis value for the whole vacuum pressure range. Adjustment Range: 1 – 50 mbar.
- Confirm your entry with **ok** to accept the change.

Menu item **Lock screen after**

You can use this menu item to specify the time after which the display is automatically locked from the last operation.

- Touch the function button to open the on-screen numeric keyboard.
- Use the keyboard to enter the desired wait time in a range of at least one (1) to a maximum of sixty (60) minutes. Entries outside this value range are not accepted!

From a continuous on-time of approx. 60 minutes, there is a risk of the display becoming burnt in and static display elements remaining visible as a visible pattern temporarily, even after switching off. It is therefore advisable to set a maximum time of 60 minutes before the screen lock is activated.

When the screen is locked, the lock icon appears on the display.

- Touch anywhere on the screen to display the screen unlock slider:

- To unlock the screen, drag the slider to the right.

Menu item **Brightness**

You can use this menu item to adjust the brightness of the display.

- Touch the function button to open the on-screen numeric keyboard.
- Define the required brightness in percent.
- Confirm your entry with **ok** to accept the change.

Menu item Max. bath temp.

In this menu, define the maximum heating bath temperature (e.g. for the processing of thermosensitive media).

- Touch the function button to open the on-screen numeric keyboard.
- Set the maximum heating bath temperature in the range between 20° C and 90° C. Entries outside this value range are not accepted!

When loading favorites and ramp profiles, it is automatically checked whether the set maximum bath temperature could affect the process. In this case, a warning message appears on the display.

- Confirm your entry with **ok** to accept the change.

Menu item Overpressure

This menu item allows you to set the value for the safety shutdown of the system when an impermissible high vacuum pressure is reached.

As soon as the vacuum pressure reaches this limit, all functions are switched off (vacuum, flask rotation, heating bath) and the flask lift is moved to the upper end position (only variants with motor lift!). In addition, a warning message appears on the display.

- Touch the function button to open the on-screen numeric keyboard.
- Define a limit value in the range between 900 mbar and 1,400 mbar (factory setting: 1,400 mbar). Entries outside this value range are not accepted!

When loading favorites and ramp profiles, it is automatically checked whether the set limit value could affect the process. In this case, a warning message appears on the display.

- Confirm your entry with **ok** to accept the change.

Menu item Leakage pressure

This menu item allows you to set the limit value for the safety shutdown of the system in case of failure to reach the leakage pressure.

This safety function is activated when the leak pressure falls below the specified minimum during the evacuation process.

If the leak pressure is not reached during evacuation (e.g. due to glass breakage or leaks), all functions are switched off (vacuum, flask rotation, heating bath) and the flask lift is moved to the upper end position (only variants with motor lift!). In addition, a warning message appears on the display.

- Touch the function button to open the on-screen numeric keyboard.
- Define a limit value in the range between 1 mbar and 1,399 mbar (factory setting: 1,399 mbar). Entries outside this value range are not accepted!
- Confirm your entry with **ok** to accept the change.

Menu item Max. pump performance

This menu item allows you to define the maximum pump performance in percent of one hundred.

By reducing the intake capacity of the pump, a higher solvent recovery rate can be achieved.

When using a speed-controlled vacuum pump, the maximum pump performance can be reduced to a maximum of 20% in 10% increments to avoid foaming or boiling delay. When using a diaphragm vacuum pump, this menu item is grayed out!

- Touch the function button to open the on-screen numeric keyboard.
- Set the maximum pump performance in the range between 20 % and 90 %. Entries outside this value range are not accepted!
- Confirm your entry with **ok** to accept the change.

Menu item Calibrate pressure sensor

The vacuum sensor is properly calibrated by the manufacturer before delivery. Normally, recalibration is not required.

PRACTICAL TIP

During a user recalibration, the evaporation flask should rotate to ensure that system does not have any leaks.

All other processes must be stopped!

- Make sure that the pressure sensor is properly connected to the device.
- Ventilate the system.
- Touch the function button to open the **Calibrate pressure sensor** mask and follow the instructions on the display:

- Use the on-screen keyboard to enter the actual atmospheric pressure.
- Close all openings on the evaporator system.
- Confirm the entry. The vacuum control starts automatically.

- Once the system is evacuated to a low level, the second calibration value can be entered: Wait until the vacuum pump does not reduce the pressure further, i.e. until the minimum possible vacuum pressure is reached.
- Enter this value and confirm the values.
- The calibration is then corrected based on the input values.

When the calibration process is canceled and the device is reset to the factory settings, the default values are loaded!

Menu item Data logging

This menu item allows you to activate/deactivate the data recording function.

- **Data logging inactive:** Touch the function button to activate the function.
- **Data logging active:** Touch the function button to deactivate the function.
- Confirm your entry with **ok** to accept the change.

For detailed information about the data recording function, see section "Data Recording" on page 120

Menu item Start process: Lift

Only devices with motor lift: In this menu, define the behavior of the flask lift at the start of the process.

- Touch the function button to switch between the two options:
 - **Start process: Lift on:** With this selection, the flask lift automatically moves to the lower end position when starting the process and the evaporation flask is thus automatically immersed in the heating bath.
 - **Start process: Lift off:** With this selection, the flask lift remains in the upper end position on starting the process; the evaporation flask must be immersed in the heating bath manually.
- Confirm your entry with **ok**.

Deactivate this function when using a foam brake and on devices using a large evaporation flask to prevent damage to the flask!

Menu item **Stop process: Heating bath**

In this menu, define the behavior of the heating bath heating at the end of the process.

- Touch the function button to switch between the two options:
 - **Stop process: Heating bath on:** With this selection, the heating bath remains switched on after the end of the process.
 - **Stop process: Heating bath off:** With this selection, the heating bath is automatically switched off at the end of the process.
- Confirm your entry with **ok**.

Menu item **Stop process: Chiller**


In this menu, define the behavior of the chiller at the end of the process.

- Touch the function button to switch between the two options:
 - **Stop process: Chiller on:** With this selection, the chiller remains switched on at the end of the process.
 - **Stop process: Chiller off:** With this selection, the chiller is automatically switched off at the end of the process.
- Confirm your entry with **ok**.

Menu item **Factory settings**

In this menu, reset the device to the factory settings (as-delivered condition). All user data (parameter adjustments, favorites, ramp profiles) will be deleted!

- Touch the **Factory settings** function button.
- The following prompt appears on the display:

- Confirm the security prompt with **ok** to restore the factory settings. The device must be restarted!
- To cancel the reset and return to the previous menu without making any changes, touch the **Cancel** button.

Selecting the user language

In the **Languages** menu (**Main menu** → **Languages**), the user language can be changed at any time. Available languages are English, German, Mandarin, French, Spanish, Portuguese, Italian, Russian and Polish.

- Touch the Languages function button.
- The list of all available menu languages opens.
- Touch the entry for the required language.
- The background color of the selected entry changes to orange.
- Press **ok** to confirm the selection and to return to the start screen of the active mode.
- To discard the entry and to return to the start window of the active mode without making any changes, touch the **Cancel** button.

Selecting formats

The following system settings can be made in the **Formats** menu (**Main menu** → **Formats**):

- Set date format (US or EU)
- Set time format (12/24-hour)
- Select heating bath fluid (water, oil)
- Set number separator (US or EU)
- Set temperature unit (degrees Celsius, Fahrenheit, Kelvin)
- Set unit for pressure display (mbar, hPa, torr)

- The active selection is highlighted in orange. Touch the desired function button to switch between options. Several changes can be made at the same time!
- Confirm the new settings with **ok**. All changes are immediately accepted without a confirmation prompt!
- To discard changes and to return to the previous menu, touch the **Cancel** button.

Data Recording

The **Data logging** function can be used to record the individual steps of a process. The function records each individual step with start and end time and stores this information in the Data logging memory.

The internal database offers a total of 30 slots with a capacity of 550,000 entries per measurement series for this.

The data is recorded at fixed intervals (ten seconds) with a fixed recording rate of 1 entry/10 seconds, i.e. the maximum recording time is approx. 152 hours.

- If the remaining capacity is 10%, a first warning message appears on the display.
- If the remaining capacity is 0%, a second warning message appears on the display and data recording is stopped. In this case, export and/or delete existing data from the system memory to record additional data.
- No running processes will be affected by the stop of the data recording!

Activating/deactivating the data logging function


The data recording function is manually activated/deactivated:

- To enable data recording, open the **Settings** menu and activate the function, see section "Menu item Data logging" on page 116.
- When the data recording function is activated, the **Start process** and **Stop process** buttons are marked with an additional icon in the footer of the start window:

Via the option **Data logging**, the recorded data can be exported to a USB stick and read out on a PC or deleted from the system.

- Open the **Main menu** and touch the function button **Data logging**.
- The Data logging list opens:

→ The recorded values show when a process was started and stopped again.

In the event of a power interruption or when the device is switched off during the current data recording, the recording process is aborted and the corresponding entry marked with "(!)". The displayed data set contains the recorded process data that was recorded up to the power interruption:

Start time: 28.09.2021 07:17	
End Time: 28.09.2021 07:17 (!)	🗑️
2 Datapoints	

Start time: 27.09.2021 16:54	🗑️
------------------------------	----

In the data log export file, the corresponding file names are also given the addition "(!)" before the file extension: "(!).csv", or "(!).json".

Deleting entries

- To delete an entry from the list, touch the associated **Delete** button. During the deletion, the assigned **Delete** button has a yellow background:


Saving entries on USB stick

Recorded entries can be saved on an USB-stick. To do this, connect a VFAT-formatted USB stick to the USB port of the display and open the data list as described.

As soon as a USB stick is connected, a **Save** button (diskette symbol) is assigned to each list entry.

- To save an entry from the list, touch the associated **Save** button. During the saving, the assigned **Save** button has a yellow background:

PRACTICAL TIP

Large amounts of data are generated during the long-term recording of process data (several days).

It is therefore recommended to save the current data recording memory on a USB stick and to remove the existing entries from the internal memory before starting the process.

For detailed information on the data recording function, see section "Data Recording" on page 120.

Connectivity

All necessary network parameters and basic settings for time synchronization can be defined in the **Hei-CONTROL Pro** menu:

- The connection status of the Hei-CONTROL Pro application is displayed in the [Status] field.

The full functionality of the Hei-CONTROL Pro application also requires the Hei-CONTROL Pro server. For more information, refer to the Hei-CONTROL Pro User Manual.

- All associated parameters can be defined manually in [static IP configuration] mode (checkbox deselected, figure on the left). Select the checkbox (figure on the right) to obtain the IP address of the device automatically.
- Select the [NTP server] checkbox to synchronize the system time of the device with a time server. This function requires a valid IP configuration!
 - Place the cursor in the corresponding input field and use the on-screen keyboard to enter the required server address.
- Select the [RS232 Distimatic Pro] checkbox if the device is to be connected directly to an automatic Hei-VOLUME Distimatic Pro module and will be controlled by it.


If the time is synchronized via a time server, the **System time/Timer** function button in the display header can only be used to select a time zone.

If an ethernet connection exists, the system time is updated automatically via the defined time server. To define the system time manually, the [NTP server] option must be deselected!

System check

All connected sensors and actuators or peripherals are displayed in the **System check** menu. Various tests can be carried out via the **System check** sub-menu to ensure the functional and operational safety of the system:

- Open the **Main menu** and touch the **System check** function button in the footer. The **System check** menu appears:

- **Fault-free** components are marked with the [OK] status indicator.
- **Faulty** components are marked with [NOK]. Please note the corresponding error messages.
- Swipe up or down on the touch panel to scroll through the list.
- Touch the **TEST** function button in the header of the **System check** menu to open the **System check** sub-menu.

CAUTION: Damage to property, loss of production

Faulty connections or signaling can lead to malfunctions or uncontrolled switching operations!

Leaks reduce the performance of the entire system.

- Correct any indicated errors and/or detected malfunctions immediately.
 - In the event of a continued fault indication/malfunction, contact your local sales representative or our technical service department, see "Contact information Heidelberg international" on page 147.
-

Cleaning the PTFE seal

The PTFE seal of the device must be cleaned after a 100-hour break-in period. To do this, follow the steps in this section.

CAUTION: Damage to the PTFE seal

Improper cleaning can damage the surfaces of the PTFE seal.

- Clean the seal surfaces with a soft, lint-free and only slightly moistened cloth.
- Never use any aggressive or abrasive cleaning agents or aids.

- Loosen the union nut and remove the condenser of the device.
- Remove the PTFE seal from the vapor tube.
- Clean the PTFE seal with a soft cloth, paying particular attention to any residue under the sealing lips.
- Put the PTFE seal back in the correct position and against the stop on the vapor tube (note the "Motor side" embossing!).
- Mount the condenser on the vapor tube again, make sure to insert the tension spring flat and do not tilt the union nut!

Operating modes

Observe the following specific safety instructions when operating the device:

WARNING

Risk of injury and/or damage to property

Due to incorrect/insufficient process settings, there is a risk that impermissible operating values are reached and, as a result, device components will be damaged and heating bath fluid and/or distillation material will escape.

- Only qualified personnel may adjust the process settings of the device.
- When determining the pressure values, observe the chemical-physical properties of the respective distillation product.
- Pay particular attention to the specific safety instructions for the processing of toxic substances!

Danger due to rotating parts

When operating the device, there is a risk of contact with rotating parts. Loose clothing, jewelry and open hair can be pulled in!

- If possible, operate the device with a guard hood (accessory) or under a closed laboratory fume hood.
- Never touch rotating parts with your hands. Do not bend over rotating parts.
- Wear tight-fitting work clothes, place jewelry before starting work, tie long hair under a hair net or under a suitable head cover.
- Use appropriate personal protective equipment (PPE).

Danger due to ejected heating bath liquid

At high rotation speeds, the bath fluid can spray away from the evaporation flask by forming film and/or swan out of the heating bath. There is a risk of scalding and/or contamination.

- If possible, operate the device with the guard hood closed.
- Adjust the speed of rotation of the evaporation flask according to the recommendations in this section.
- Use appropriate personal protective equipment (PPE).
- Clean dirty surfaces immediately!


After switching on the device and initializing the software, the start page of the last active operating mode appears.

Open the **Main** menu and touch the **Applications** function button (see section "Applications" on page 111). Activate the **Base** mode (see section "Base mode" on page 127) or the **Advanced** mode (see section "Advanced mode" on page 128).

Base mode

After switching on the device and initializing the software, the start page of the last active operating mode appears. Open the **Main menu** and touch the **Applications** function button (see section "Applications" on page 111). Activate the **Base mode**.

In **Base mode**, the operating parameters vacuum pressure, rotational speed and heating bath temperature are entered manually by the operator:

- Touch a parameter's field to open the on-screen keyboard and set a desired value. The input values are accepted without a security prompt:
 - Coolant temperature (setting range: -10 – 20 °C)
 - Vacuum pressure (setting range: 1 – 1,400 mbar)
 - Rotation speed (10 – 280 rpm)
 - Heating bath temperature (20 – 80 °C)
- Touch the icons of the parameter displays **coolant temperature, vacuum pressure, rotation** and **heating bath temperature** to activate (icon changes from white to orange) or to deactivate (icon changes from orange to white) the assigned process independently of the other processes.
- Touch the **Start Process** button to start all processes with the displayed setpoints at the same time.
- Touch the **Stop Process** button to stop all ongoing processes at the same time.

Advanced mode

After switching on the device and initializing the software, the start page of the last active operating mode appears. Open the **Main menu** and touch the **Applications** function button (see section "Applications" on page 111). Activate the **Advanced** mode.

The following modes are available in **Advanced** mode:

- Dynamic AUTOaccurate
- Ramp profiles
- Favorites mode

Select the desired option via the assigned **Mode** icon in the header of the start page (in the example, Dynamic AUTOaccurate):

Dynamic AUTOaccurate (DAA) mode

In Dynamic AUTOaccurate (DAA) mode, the vacuum pressure is controlled depending on the temperature.

REQUIREMENTS

- For operation in Dynamic AUTOaccurate (DAA) mode, the device has to be equipped with an optional AUTOaccurate sensor.
- The AUTOaccurate sensor can only be used in combination with glassware G3 or G6!

The temperature-dependent vacuum control in DAA mode is based on two values:

- **Condenser temperature:** Temperature reading on the AUTOaccurate sensor in the condenser
- **Intensity:** Relative temperature rise to the closing of the vacuum valve, based on the measured condenser temperature T(AA)

At the start of the process, the condenser temperature T(AA) is measured using the AUTOaccurate sensor. When reaching a stable level, the actual evaporation process starts. This "pre-process" can be manually canceled at any time. In this case, the evaporation process is started immediately with the currently measured condenser temperature value.


During the evaporation process, the temperature at the AUTOaccurate sensor increases relative to the initial condenser temperature T(AA). If this temperature rise reaches the defined intensity (parameter **Intensity**), the evaporation process is stopped.

As soon as the temperature at the AUTOaccurate sensor has dropped by a defined value, the system is evacuated and the evaporation process starts again. If the temperature rise at the AUTOaccurate sensor reaches the defined intensity again (at relatively higher vacuum pressure), the vacuum pump is switched off again, etc.

In this way, a consistent performance is achieved in the processing of solvents or mixtures throughout the distillation process.

During multiple solvent distillation, the described process steps are repeated until a defined maximum vacuum pressure is reached (input value **End pressure**). The distillation process is then switched off.

→ Activate the **Advanced** mode and change to **Dynamic AUTOaccurate**:

Touch a parameter's field to open the on-screen keyboard and set a desired value. The input values are accepted without a security prompt:

- Coolant temperature (setting range: -10 – 20 °C)
- Rotation speed (10 – 280 rpm)
- Heating bath temperature (20 – 80 °C)

→ Switch to the second page of the screen (swipe left on the touch screen) and set the type of distillation process, the intensity and a final pressure (for multiple solvent distillation only):

The value **Switch point Ts** results from the measured condenser temperature $T(AA)$ and the defined intensity.

This value is displayed on the second screen after the process is started, see figure.

If the intensity is adjusted during the ongoing DAA process, the value must be changed by at least 0.5 °C. Otherwise, the change will not be taken into account.

- Go back to the home page. The following options are available:
 - Touch the **coolant temperature, rotation and heating bath temperature** icons to activate (icon changes from white to orange) or to deactivate (icon changes from orange to white) the assigned process independently of the other processes.
 - Touch the **Warm Up** button to activate the heating function for the heating bath.
 - Touch the **Start Process** button to start all processes with the displayed setpoints at the same time.

Use the **Start process: Lift on/off** parameter to define the behavior of the motor lift at the start of the process, see section "Menu item Start process: Lift" on page 116.

- Touch the **Stop Process** button to stop all ongoing processes at the same time.

Use the **Stop process: Heating bath on/off** parameter to define the behavior of the heating bath heater on stopping the process: Heating bath on = heating bath heating is not switched off when stopping the process! See section "Menu item Stop process: Heating bath" on page 117).


Use the **Stop process: Condenser on/off** parameter to define the behavior of the condenser on stopping the process: Chiller on = chiller is not switched off when the process is stopped! See section "Menu item Stop process: Chiller" on page 117.

Ramp profiles mode

In Ramp profiles mode, user-specific process sequences can be stored as a ramp profile. Activate the **Advanced** mode and change to **Ramp profiles**:

- The last activated ramp profile is loaded automatically. Switch from the home screen to the second screen (swipe left on the touch screen). The following options are available:
- Load ramp profile (**Upload** button)
 - Create/edit ramp profile (**New** button or Upload option)
 - Delete ramp profile (**Delete** button)
 - Save ramp profiles on USB stick (**List to USB** button – available as soon as a USB stick is inserted on the control panel!)

Load ramp profile

- If you have already saved a profile, highlight an existing profile in the list displayed. The list can contain several screens: Use the two arrow keys to move between the list pages.
- Press **Upload** to load the profile and to return to the home page.


The designation of the loaded ramp profile is displayed in the process visualization area.

Screen page 4 in **Ramp profiles** mode shows the associated diagrams for vacuum pressure, rotation speed, and condenser temperature. The graphs are touch-sensitive: Touch a chart to zoom in and out.

- Touch the **Start process** button to process the ramp profile.

Create ramp profile

- Touch the **New** button to create a new rap profile:

- Use the **[+]** key to add a new line:
- Touch the **hh:mm** field in the new line and use the on-screen keyboard to specify the duration of the process step (setting range 00:01 to 99:9).
- Touch the **mbar** field in the new line and define the vacuum pressure for the process step (setting range 1 – 1,400 mbar).
- Touch the **rpm** field in the new line and set the rotation speed for the process step (setting range 10 – 280 rpm).
- Touch the **°C** field in the new line and set the heating bath temperature for the process step (setting range 20 – 210 °C).

For each individual value in the list (vacuum pressure [**mbar**], rotation speed [**rpm**] and heating bath temperature [**°C**]), it can be individually determined whether the setpoint should be reached abruptly or gradually linearly.

Touch the symbol assigned to toggle between the two options:

Option abruptly:

Option gradually linear:

- Use the **[+]** key to add additional lines if necessary and set the parameters for each process step as described.
- To delete a line, highlight it and touch the **[-]** button.
- After entering all required lines and values, open the on-screen keyboard by touching the **Save** button.
- Assign an individual designation to the ramp profile.

It is not possible to assign designations for ramp profiles twice in the system. If a desired designation already exists, a warning message appears on the display.

In this case, change the designation or overwrite the existing ramp profile. This operation cannot be undone!

- Confirm the new designation to store the ramp profile in the system.
- Switch to screen 5 in **ramp profiles** mode to define the individual process-end-parameters. Options:
 - **Stop vacuum / Maintain vacuum** (touch the entry to switch between the two options)
 - **Hysteresis** (touch the entry to open the automatic/manual hysteresis selection menu, see also section "Menu item Hysteresis Δp " on page 113)
 - **Stop rotation / Maintain rotation** (touch the entry to switch between the two options)


- **Stop heating bath / Hold heating bath** (touch the entry to switch between the two options)
- **Boot lift / Hold lift** (touch the entry to switch between the two options)
- **Stop chiller / Maintain chiller** (touch the entry to switch between the two options)

The settings made here are prioritized and suppress the general parameter settings in the **Settings** menu!

Edit ramp profiles

- Load the desired ramp profile as described above.
- Switch to the parameters page (screen 3 in **ramp profiles** mode): Swipe left twice from the start page of the loaded ramp profile on the touch panel:

- From the list of all process steps of the ramp profile, select all the affected lines one after the other (highlighted in orange) and edit the values as described above (see section "Touch the Start process button to process the ramp profile." on page 132).
 - You can change several process steps at the same time.
 - The designation of a modified ramp profile is marked with an asterisk (*).
- If necessary, adjust the end-of-process settings on screen 5 in the **ramp profiles** mode, see section "Create ramp profile" on page 132.
- Confirm all changes with **Save**. The on-screen keyboard opens.
- Confirm the designation to store the ramp profile with the changed values in the system.

Touch the **Cancel** button to cancel the operation at any time and at any level. All changes made to the ramp profile before are discarded.

- Confirm the additional confirmation prompt with **ok** to accept all changes.

Alternatively, assign an individual designation to the modified ramp profile to create a new one with similar values based on the existing one, see also section "Touch the Start process button to process the ramp profile." on page 132.

Delete ramp profiles

- Open the list of ramp profiles (see section "Load ramp profile" on page 132).
- Highlight an existing profile. The list can contain several screens: Use the two arrow keys to move between the list pages.
- Touch the **Delete** button and confirm the security prompt to permanently delete the profile.

Assign quick access buttons

- Load a desired profile (see section "Load ramp profile" on page 132).
- In the footer of the home page, touch one of the two quick access buttons for about two seconds:

- As soon as the icon color changes from white to orange, the quick access buttons is assigned to the previously loaded ramp profile (in the example, quick access button 1).

Export ramp profiles

- As soon as a USB stick is inserted into the control panel, the **List to USB** button is displayed in the footer of the ramp profiles view.
- Touch the **List on USB** button to save all listed ramp profiles on the connected USB stick.

Import ramp profiles

- As soon as a USB stick with stored ramp profiles is inserted into the control panel, a dialog box appears prompting you to import the existing data.
- To import the existing data, confirm with **ok**. Alternatively, close the dialog box without importing by touching the **Cancel** button.


Only ramp profiles can be imported in the **ramp profiles** mode. Unrecognized formats (e.g. favorites) are not displayed.


Files with file names that already exist in system memory cannot be imported. In such a case, you can rename the file stored in the system before the import or overwrite this file with the import file.

Favorites mode

In favorites mode, user-specific parameter settings for the distillation of certain substances can be stored as favorites. Activate the **Advanced** mode and change to **favorites**:

- The last activated parameter set is loaded automatically. Switch from the home screen to the second screen (swipe left on the touch screen). The following options are available:
- Load favorite (**Upload** button)
 - Create/edit favorites (**New** button or Upload option)
 - Delete favorites (**Delete** button)
 - Save favorites on USB stick (**List on USB** button – available as soon as a USB stick is inserted on the control panel!)

Load favorite


- If you have already saved favorites, highlight an existing favorite (process parameter record) in the list displayed. The list can contain several screens: Use the two arrow keys to move between the list pages.
- Press **Upload** to load the favorite and to return to the home page.
- Touch the **Start process** button to process the favorite.

The process must be stopped manually or by timer in this mode. There is no automatic switch-off!

Create favorites

- Touch the **New** button to create a new record:

- Switch to screen 3 in **favorites** mode to define the individual parameters. Options:
 - **Vacuum** (vacuum pressure, setting range 1 – 1,400 mbar)
 - **Hysteresis** (touch the entry to open the automatic/manual hysteresis selection menu, see also section "Menu item Hysteresis $\Delta\rho$ " on page 113)
 - **Rotation** (rotation speed, setting range 10 – 280 rpm)
 - **Heating bath** (heating bath temperature, setting range 20 – 210 °C)
 - **Chiller** (chiller temperature, setting range -10 – 20 °C)
- After entering all required values, press **Save** to open the on-screen keyboard
- Assign an individual designation to the favorite.

It is not possible to assign designations for favorites twice in the system. If a desired designation already exists, a warning message appears on the display.

In this case, change the designation or overwrite the existing favorite. This operation cannot be undone!

Edit favorites

- Load the desired favorite as described above.
- Switch to the parameters page (screen 3 in **favorites** mode): Swipe left twice from the start page of the loaded favorite on the touch panel:
- To edit the individual parameters, follow the steps in the previous section.
- Confirm all changes with **Save**. The on-screen keyboard opens.
- Confirm the designation to store the favorite with the changed values in the system.

Touch the **Cancel** button to cancel the operation at any time and at any level. All changes made to the ramp profile before are discarded.

- Confirm the additional confirmation prompt with **ok** to accept all changes.


Alternatively, assign an individual designation to the modified favorite to create a new one with similar values based on the existing one.

Delete favorites

- Open the list of favorites (see section "Load favorite" on page 138).
- Highlight an existing favorite. The list can contain several screens: Use the two arrow keys to move between the list pages.
- Touch the **Delete** button and confirm the security prompt to permanently delete the favorite.

Assign quick access buttons

- Load a desired favorite (see section "Load favorite" on page 138).
- In the footer of the home page, touch one of the two quick access buttons for about two seconds:

- As soon as the icon color changes from white to orange, the quick access buttons is assigned to the previously loaded favorite (in the example, quick access button 1).

Export favorites

- As soon as a USB stick is inserted into the control panel, the **List to USB** button is displayed in the footer of the ramp profiles view.
- Touch the **List to USB** button to save all listed favorites to the attached USB stick.

Import favorites

- As soon as a USB stick with stored favorites is inserted into the control panel, a dialog box appears prompting you to import the existing data.
- To import the existing data, confirm with **ok**. Alternatively, close the dialog box without importing by touching the **Cancel** button.

Only favorites can be imported in the **favorites** mode. Unrecognized formats (e.g. ramp profiles) are not displayed. Files with file names that already exist in system memory cannot be imported. In such a case, you can rename the file stored in the system before the import or overwrite this file with the import file.

Troubleshooting

Failure	Possible cause – remedy
Device does not turn on Display remains dark after switching on	No power supply, check power supply cord.
	Fuse defective, replace.
	On/off switch to OFF, switch on.
Device does not heat up	On/off switch defective, contact technical service.
	Connection cable heating bath is not connected, connect.
	Heating bath defective, contact technical service.
Drive does not turn	Overheat protection has tripped: Let the heating bath cool down and empty, reset the overheating protection.
	No rotation movement, press speed controller.
	Drive defective, contact technical service.
Engine lift not running	Lift in end position, move in opposite direction.
	Lift at height stop, adjust height stop.
	Mechanical defect/motor defective, contact technical service.
No vacuum	Vacuum pump on/off switch OFF, switch on
	Vacuum valve defective, replace
Vacuum insufficient	System leaking, check seals and connections, check joints, grease if necessary.
	Vacuum pump defective, observe manufacturer's instructions for vacuum pump.
Device shuts down unexpectedly	Fuses defective, replace.
	Emergency stop when detecting pressure (see settings, overpressure and leakage pressure), check system and switch on the device again, calibrate the vacuum sensor.
Boiling temperature is not displayed	Boiling temperature sensor not connected, connect or check connection, switch device off and on again if necessary.
	Boiling temperature sensor defective, visual inspection for defect, contact technical service, replace sensor if necessary.
	Temperature outside the measuring range, heat the sensor or let it cool down.

Automatic mode (DAA) cannot be selected	<p>AUTOaccurate sensor not connected, connect or check the connection, switch off and on the device if necessary.</p> <p>AUTOaccurate sensor defective, visual inspection for defect, contact technical service, replace sensor if necessary.</p> <p>Temperature outside the measuring range, heat the sensor or let it cool down.</p>
Control box not connected or defective	Check connection between control box and evaporator, contact technical service.
Ventilation valve defective	Short circuit in the ventilation valve, contact technical service.
Vacuum valve not connected or defective	<p>Vacuum valve not connected or connection between control box and vacuum valve interrupted, connect vacuum valve, check connection between control box and vacuum valve, contact technical service.</p> <p>Short circuit in the vacuum valve, contact technical service.</p>
Vacuum pump not connected or defective	Vacuum pump not connected or connection between Control Box and vacuum pump interrupted, connect vacuum pump, check connection between Control Box and vacuum pump.

All system messages (error messages, warnings, information) are shown in plain text on the display.

Follow the instructions on the display.

In case of recurring errors, please contact the responsible sales department or our technical service. Contact address see "Contact information Heidolph international" on page 147.

Technical specifications

General device data	
Model	Rotary evaporator Hei-VAP Expert/Ulimate
Dimensions (W × H × D)	Variant with motor driven lift & glassware G3: 739 × 887 × 477 mm
	Variant with manual lift & glassware G3: 739 × 887 × 532 mm
Weight	approx. 15 kg, without glassware
Acoustic pressure	< 85 (dB(A)) (in accordance with IEC 61010)
Flask drive	EC motor
Speed range	10 – 280 rpm
Stroke flask lift	155 mm
Lifting speed flask lift	55 mm/s
Evaporation rates (L/h) ΔT^* 40 °C (in continuous operation)	Toluene 8.5
	Acetone 5.8
	Ethanol 3.5
	Water 1.2
Maximum cooling surface	1.4 m ²
Protection class	Device IP20
	Control panel IP42
	Control-Box IP42
	Connection cable heating bath IP67
Electrical data	
Rated voltage	230 V (50/60 Hz) or 115 V (50/60 Hz)
Connection	L+N+PE
Protection class	I  (IEC 61140)
Overvoltage category	II
Degree of pollution	2
Power input	1,400 W max.

Heating bath, cooling	
Heating bath	Inner Ø: 253 mm, outer Ø: 291 mm Material: V4A steel (1.4404) X2CrNiMo 17-12-2 Volume: 4.5 L
Heating power	1,300 W
Temperature range heating bath	20 – 100 °C H ₂ O / 20 – 210 °C oil
Overtemperature protection	Cut-off at 5 °C deviation from the set temperature
Overheat protection	Threshold value 250 °C
Control	electronic
Control accuracy	±1 °C
Permissible ambient conditions	
Operating temperature	5 °C – 31 °C at up to 80 % rel. humidity 32 °C – 40 °C at up to 50 % rel. humidity (decreasing linearly)
Installation altitude	up to 2,000 m asl

Scope of delivery

Item	Quantity	Product no.
Rotary evaporator Hei-VAP Expert, complete*	1	see package list
Rotary evaporator Hei-VAP Ultimate, complete*		
Operating instructions	1	01-005-005-86
Warranty registration	1	01-006-002-78

***The actual delivery depends on the respective device variant.**

Accessories

Detailed information on the available accessories for your device variant can be found on our website at www.heidolph-instruments.com.

In case of need, contact an authorized dealer or our technical service, see "Contact information Heidolph international" on page 147.

Device service

When carrying out service work on the device (cleaning, maintenance, repair), observe the general instructions and safety information described in this section.

WARNING: Danger of electric shock

Live components are installed inside the device.

When opening the device, there is a risk of touching live components.

- Switch the device's main switch off and disconnect it from the power supply before carrying out maintenance work, cleaning, or repairs.

Penetrating liquid poses the danger of an electric shock.

- When cleaning, avoid the penetration of liquids.

General cleaning instructions

Wipe all surfaces and the control panel with a damp cloth if necessary. Persistent contamination can be removed with mild soapy water.

CAUTION: Damage to the device

Improper cleaning can damage the surfaces of the device.

Penetrating liquid can damage the electronic components inside the device.

- Clean the device's surfaces with a soft, lint-free and only slightly moistened cloth.
- Never use any aggressive or abrasive cleaning agents or aids.

Repairs

Repairs to the device may only be carried out by authorized experts! Unauthorized repairs during the warranty period will result in the loss of the warranty claim.

The owner is solely liable for damage caused by unauthorized repairs.

In case of a need for repair, contact an authorized dealer or our technical service, see "Contact information Heidolph international" on page 147.

Include the completed certificate of decontamination with every device return, see "Certificate of decontamination" on page 148.

Maintenance

There are no user-serviceable components in the unit housing. If necessary (in the event of abnormal operating behavior such as excessive noise or heat generation, for example), contact our technical service, see "Contact information Heidolph international" on page 147.

- Check for tight fit the caps and fasteners of all glass components at regular intervals and especially before each start-up of the device.
- Check the PTFE seal of the EASY LOCK flask clamping system at regular intervals for heavy contamination and visible damage.
- Clean the PTFE seal of the EASY LOCK flask clamping system in regular preventive manner to prevent excessive deposits or encrustations.

CAUTION: Risk of property damage, loss of production

At loose/leaky connections, ambient air is sucked in during the vacuum build-up and the required vacuum pressure cannot be reached in a stable manner!

In case of excess pressure, there is a risk that fluid will escape from the system! The device and surrounding areas can be severely contaminated by leaking media!

- Tighten loose caps/fasteners by hand or correct the seating of the retaining clips.
 - Clean the PTFE seal of the EASY LOCK flask clamping system at regular intervals and depending on the use of the device and the processed media!
 - Replace the PTFE seal of the EASY LOCK flask clamping system if it shows visible damage!
-

Dismantling

Observe the instructions given in the associated assembly instructions when dismantling the device.

Disposal

- When disposing of the device, observe the provisions of the WEEE Directive 2012/19/EU and its transposition into national law in the country of use.
 - When disposing of portable batteries, observe the provisions of the European Battery Directive 2013/56/EU and their implementation in national law in the country of use.
 - Check the device and all components for residues of substances that are hazardous to health, the environment and biohazardous before disposal.
 - Properly remove and dispose of residues of substances that are hazardous to health, the environment and are biohazardous!
-

Contact information Heidolph international

Heidolph Instruments North America

Phone: 1-866-650-9604
E-mail: service@heidolph.com
www.heidolphNA.com

Heidolph Instruments United Kingdom

Phone: 01799 - 5133-20
E-mail: service@radleys.co.uk
www.heidolph-instruments.co.uk

Local distributors

To find your local distributor please visit www.heidolph.com

Warranty Statement

Heidolph Instruments provides a three-year warranty against material and manufacturing defects.

Glass and wear parts, transportation damage, and damage resulting from improper handling or non-intended use of the product are excluded from the warranty.

The warranty period for registered products begins on the date of purchase. Register the product with the enclosed warranty card or on our homepage www.heidolph.com.

For non-registered products, the warranty period begins with the date of the serial production (to be determined by the serial number).

In the event of material or manufacturing defects, the product will either be repaired or replaced free of charge within the warranty period.

Certificate of decontamination

Enclose the declaration of no objection, duly completed, with your device return.
Submissions without a declaration of no objection cannot be processed!

Certificate of decontamination

IN CASE OF RETURNS

Please fill in the required fields.

Note: The sender must package the goods properly and appropriately for transport.

Heidolph Instruments GmbH & Co. KG
Walpersdorfer Straße 12
91126 Schwabach, Germany
Phone: +49 (0) 9122 9920-380
Fax: +49 (0) 9122 9920-19
Email: service@heidolph.de

SENDER

Name _____ First name _____
Company/institution _____ Department _____
Workgroup _____
Address _____
ZC/City _____
Country _____ Phone _____
Email _____

DEVICE DETAILS

Article number _____ Serial no. _____
Ticket number _____
Reason for sending in _____

Has the device been cleaned, decontaminated/disinfected? Yes No (Please mark as applicable)
If yes, which measures were carried out?

Does this device pose a risk to people and/or the environment due to the processing of substances that are hazardous to health, the environment and/or are biohazardous? Yes No (Please mark as applicable)
If yes, with which substances did the device come into contact?

LEGALLY BINDING DECLARATION

The principal/consignor is aware that they are liable to the agent/consignee for losses or damage incurred due to incomplete and incorrect information.

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hereby declare, that the product of the brand "Heidolph" (from serial number on 200067818) designated below is in compliance with the basic requirements of all applicable EU-directives stated below with regard to design, type of model sold and manufactured by us. This certificate will be invalid if the product is modified without the prior written consent and agreement of the manufacturer.

Hei-VAP Core	517-61000-00, 571-00000-00, 572-00000-00 569-50000-00, 569-51000-00
Hei-VAP Expert	517-61000-00, 573-00000-00, 574-00000-00 569-60000-00, 569-61000-00
Hei-VAP Ultimate	517-61000-00, 573-00000-00, 574-00000-00 569-70000-00, 569-71000-00
Hei VAP Expert Control	517-61000-00, 573-00000-00, 574-00000-00 569-60000-00, 569-61000-00, 569-00200-00
Hei-VAP Ultimate Control	517-61000-00, 573-00000-00, 574-00000-00 569-70000-00, 569-71000-00, 569-00200-00

Maschinenrichtlinie / Machinery Directive 2006/42/EG
EMV-Richtlinie / Electromagnetic Compatibility Directive 2014/30/EU
Delegierte (RoHS-) Richtlinie / Delegated (RoHS) Directive 2015/863/EU
Angewandte (harmonisierte) Normen / (Harmonized) Standards applied:
EN ISO 12100:2010, EN 61326-1:2013, EN 61010-1:2010, EN IEC 63000:2018

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UK CA

Declaration of Conformity

In accordance with UK Government guidance

This declaration is issued under the sole responsibility of the manufacturer,
Heidolph Instruments GmbH & Co. KG
Walpersdorfer Straße 12
91126 Schwabach / Germany

Product: Laboratory benchtop evaporator

Model:

Hei-VAP Core	517-61000-00, 571-00000-00, 572-00000-00, 569-50000-00, 569-51000-00
Hei-VAP Expert	517-61000-00, 573-00000-00, 574-00000-00, 569-60000-00, 569-61000-00
Hei-VAP Ultimate	517-61000-00, 573-00000-00, 574-00000-00, 569-70000-00, 569-71000-00
Hei VAP Expert Control	517-61000-00, 573-00000-00, 574-00000-00, 569-60000-00, 569-61000-00, 569-00200-00
Hei-VAP Ultimate Control	517-61000-00, 573-00000-00, 574-00000-00, 569-70000-00, 569-71000-00, 569-00200-00

Description:

Rotary Evaporators with heating bath and human machine interface

The object of the declaration described above is in conformity with the relevant
UK Statutory Instruments (and their amendments):

2008 No. 1597

The Supply of Machinery (Safety) Regulations

2008 2016 No. 1091

The Electromagnetic Compatibility Regulations 2016

2012 No. 3032

*The Restriction of the Use of Certain Hazardous Substances in
Electrical and Electronic Equipment Regulations 2012*

and complies with the following technical standards :

EN ISO 12100:2010, EN 61326-1:2013, EN 61010-1:2010, EN IEC 63000:2018

UK Authorised Representative (for authorities only):

ProductIP (UK) Ltd.
8. Northumberland Av.
London WC2N 5BY

Signed for and on behalf of Heidolph Instruments GmbH & Co. KG
Walpersdorfer Straße 12, 91126 Schwabach / Germany

Schwabach, 01.09.2021

Wolfgang Jaenicke
Managing Director

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With this declaration, we confirm (according to current knowledge) that all sold laboratory devices by Heidolph Instruments GmbH & Co. KG fulfill the requirements of the EU directive 2011 / 65 / EU (RoHS) and the amended of directive 2015 / 863.

All devices are compatible with the requirement of the RoHS for the following materials:

Max. 0,01% of the weight in homogeneous material for cadmium and max. 0,1 % of the weight in homogeneous material for lead, mercury, hexavalent chromium, polybrominated biphenyl, polybrominated diphenyl ether, Di (2-ethylhexyl) phthalate, butyl benzyl phthalate, dibutyl phthalate, diisobutyl phthalate.

In the case of individual assemblies, maximum concentrations maybe exceeded within the permissible exceptions of the Directive.

Schwabach, 22.07.2019

Stefan Peters
Vice President Marketing, Innovation & Technologie

Marcell Sarré
Vice President Quality Management & Technical Service

China RoHS Certification

China RoHS DECLARATION OF CONFORMITY

Heidolph Instruments GmbH & Co.KG has made reasonable efforts to ensure that hazardous materials and substances may not be used in its products.

In order to determine the concentration of hazardous substances in all homogeneous materials of the subassemblies, a "Product Conformity Assessment" (PCA) procedure was performed. As defined in GB/T 26572 the "Maximum Concentration Value" limits (MCV) apply to these restricted substances:

- Lead (Pb): 0.1%
- Mercury (Hg): 0.1%
- Cadmium (Cd): 0.01%
- Hexavalent chromium (Cr(VI)): 0.1%
- Polybrominated biphenyls (PBB): 0.1%
- Polybrominated diphenyl ether (PBDE): 0.1%

Environmental Friendly Use Period (EFUP)

EFUP defines the period in years during which the hazardous substances contained in electrical and electronic products will not leak or mutate under normal operating conditions. During normal use by the user such electrical and electronic products will not result in serious environmental pollution, cause serious bodily injury or damage to the user's assets.

The Environmental Friendly Use Period for Heidolph Instruments GmbH & Co.KG products is 25 years.

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This table is created according to SJ/T 11364-2014

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Note: Table applies to all products. Some of the components or parts listed above may not be part of the enclosed product.

- O: 表示该有毒有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下。
- O: Indicates that the above mentioned hazardous substance contained in all homogeneous materials of the part is below the required limit as defined in GB/T 26572.

- X: 表示该有毒有害物质至少在该部件某一均质材料中的含量超出GB/T 26572规定的限量要求。
- X: Indicates that the above mentioned hazardous substance contained in at least one of the homogeneous materials of this part is above the required limit as defined in GB/T 26572.

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Apart from the disclosures in the above table, the subassemblies are not intentionally manufactured or formulated with lead (Pb), mercury (Hg), cadmium (Cd), hexavalent chromium (CrVI), polybrominated biphenyls (PBB), and polybrominated diphenyl ethers (PBDE).

Products manufactured by Heidolph Instruments GmbH & Co.KG may enter into further devices or can be used together with other appliances .

With these products and appliances in particular, Heidolph Instruments GmbH & Co.KG will not take responsibility for the EFUP of those products and appliances.

Schwabach, 06.08.2021

Wolfgang Jaenicke
Chief Executive Officer CEO

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Subject to change without notice. The printed version of this document is not regularly updated. The latest issue of this document can be found by visiting our homepage.

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MATERIAL CONTENT DECLARATION FOR Heidolph Instruments GmbH & Co. KG PRODUCTS
© Heidolph Instruments GmbH & Co. KG