

## **Distillation**

Distillation is a process by which a mixture of substances are separated thermally. It utilises the evaporation of the liquid sample by heating, and then condensation of its vapours by cooling and subsequently collected in a receiver.

## **Theoretical background**

Distillation (Latin; destillare, to drip down) is the evaporation of a liquid consisting of any number of components and the subsequent condensation of the vapour, to give the distillate. Fractionated distillation is the thermal separation of at least two distillates on the basis of them having different boiling points.

## Method

Materials which are temperature sensitive must be treated with care to avoid decomposition. Therefore, is it important to have gentle, indirect heating. The working principle of the BUCHI Glass Oven B-585 utilises a glass tube which is coated with a transparent and electrically conductive layer. This way of heating delivers moderate heat radiation, which is evenly distributed across the entire heating surface and facilitates heat entry of the material without direct contact of the heat source, thus minimising thermal decomposition. The glass on which the layer is coated on has good heat conductivity and high thermal capacity, therefore the majority of the heat generated is radiated inward.

The electrically conductive layer provides fast heating to the oven and since 75% of the heat is supplied to the material is by radiation and only 25% by convection, sample heating is quick. The oven heats up to 300°C and the internal temperature is monitored by an integrated sensor. And since the glass oven is completely transparent, visual observation of the sample is possible, thereby any decomposition can be observed and actions can be taken.

Two type of distillation are possible with the glass oven: simple distillation and fractionated distillation. Simple distillation is accomplished by having one ball in the oven and another ball outside, being cooled. Fractionated distillation is starts with up to balls in the oven, and one outside. As each fraction is collected, one of the balls in the oven is shifted outside to collect the next fraction. The process of simple distillation can be compared to that of a rotary evaporator, but on a much smaller scale.



## B-585

# Glass Oven



## Simple distillation with the B-585 ball tube

The distillation is achieved by placing the sample in a glass ball, positioned inside the oven. A second ball, which is connected to the first one, is placed outside the oven and is cooled using a cooling device. The oven is then heated to enable distillation. The evaporated solvent or product is condensed by the cooling device, preventing it from being sucked into the pump and subsequently, the atmosphere. Rotation of the ball tube optimises heat transfer and reduces boiling of sample.

The heating temperature, which is set at around 10°C to 40°C above the boiling temperature of the sample (depending upon type of material and pressure), is regulated by a vacuum controller, are the relevant. The heating temperature is typically set at around 10°C to 40°C above the boiling temperature of the sample (depending on type of material and pressure). A vacuum controller may be used to regulate process parameters for distillation.

For thermally unstable substances, it is preferable to use low vacuum and for solvents with very low boiling point, it is advisable to use a vacuum system with a secondary condenser in order to condense residual sol-

## Fractionated distillation


Mixtures of products can be separated as long as there is a difference of 20°C or more between the individual components. As many balls as components in the mixture can be interconnected (but only a maximum of four balls may be placed in the oven.

Example: A three-component mixture needs 1 end-ball and 2 mid-balls. The mixture that is to be separated is placed into the end ball. All the balls except for the last empty glass ball are inserted into the oven chamber, with the last glass ball cooled in the cooling device (with ice or dry ice-alcohol mixture).

The temperature of the oven then is increased until distillation of the component with lowest boiling point begins.

Once the distillation of this fraction is complete, i.e. volume does not increase in the received fraction (in the cooled ball), the distillation is halted.

The next ball is then pulled from the oven and the temperature is increased. Thus the next high-boiling component is distilled and deposited in the second cooled ball. This procedure repeats itself until all balls are filled with fraction.





Quality in your hands

**BÜCHI** Labortechnik AG Postfach 9230 Flawil Schweiz Tel. +41 71 394 63 63 Fax +41 71 394 65 65 buchi@buchi.com www.buchi.com

**BUCHI** Analytical Inc. 19 Lukens Drive New Castle Delaware 19720 USA Tel. +1 302 652 3000 Fax +1 302 652 8777

us-sales@buchi.com

www.buchi-analytical.com

**BÜCHI** Labortechnik GmbH Postfach 10 03 51 45003 Essen Deutschland Freecall 0800 414 0414 Fax +49 201 237 082 deutschland@buchi.com www.buechigmbh.de

BUCHI Analytical Ltd. Link House Stakehill Industrial Park Middleton Manchester M24 2RH Great Britain Tel +44 161 654 60 60 Fax +44 161 643 75 10 uk-sales@buchi.com www.buchi-analytical.co.uk **BÜCHI** Labortechnik GmbH Branch Office Netherlands Postbus 142 3340 AC Hendrik Ido Ambacht The Netherlands Tel. +31 78 684 94 29 Fax +31 78 684 94 30 netherlands@buchi.com www.buchi.com

Nihon BUCHI K.K. 7F Kojima Bldg. 2-11-10 Ueno Taito-ku Tokyo 110-0005 Japan Tel. +81 3 5807 5599 Fax +81 3 5807 5598 nihon@buchi.com www.buchi.com

BÜCHI Italia s.r.l Centro Direzionale, Milano Fiori Pal A-4 Strada 4 20090 Assago (MI) Italia Tel. +39 02 824 50 11

Fax +39 02 57 51 28 55 italia@buchi.com www.buchi.it

**BUCHI** Sarl 5. rue du Pont des Halles 7 A du Delta 94656 Rungis Cedex France Tél. +33 1 56 70 62 50 Fax +33 1 46 86 00 31 france@buchi.com

www.buchi.fr

**BUCHI** SMP Services Private Ltd. 201, Magnum Opus Shantinagar Industrial Area Vakola, Santacruz (East) Mumbai 400 055 India Tel. +91 22 569 89 450 Fax +91 22 569 89 452

smplisp@vsnl.com

www.buchi.com

BACC BUCHI

Thailand

**ASEAN Competence Center** 

c/o Becthai Bangkok

300 Phaholyothin Road

Samsennai, Phayathai Bangkok 10400 Tel. +66 1 829 49 10 bacc@buchi.com www.buchi.com

BÜCHI Labortechnik AG 94152 0304 2000 Technical subject to alterations. Quality System ISO 9001